

Fig. S1




Fig. S2




Fig. S3

## Document S1

### System requirements

A Mac with a 64-bit processor is required. Windows and Linux are unsupported. Recommended and tested systems are as follows.

	Recommended	Tested
OS	Mac OS 10.8 or later	10.9.1, 10.8.4, 10.7.5, 10.6.8 <sup>a</sup>
CPU	Intel Core i5 or later	Intel Core i5, Core 2 Duo
RAM	4GB or more	8, 4, 2 <sup>b</sup> GB

a, Maybe incompatible; b, too slow.

Prior to establishing the initial settings of SraTailor, Xcode (Command line tools) and Aspera Connect must be installed; both are freely available. Installation of these software packages is described in the SraTailor instruction manual, which is available on our web site (<http://www.dev.med.kyushu-u.ac.jp/manual/>).

### Supported dataset types

SraTailor can process standard sequence data that are directly mapped onto a specified genome assembly. Therefore, SraTailor does not accept SRAs derived from reads containing technical sequences (e.g. barcode and index tags) or intra-genomic ligations (e.g., chromosome conformation capture experiments). However, technical sequences may be removed by setting Bowtie options (-5 or -3). SraTailor automatically discerns single or paired reads for a specified GSM number by accessing the corresponding metadata web page. Paired reads are mapped by bowtie2 in paired-end mode (-1 and -2 inputs) to produce the SAM and BAM files; however, in the resultant BigWig, BedGraph, and peak files, the inserts of the paired ends are excluded. Colorspace sequences used in the SOLiD platform are unsupported.